

Applications Actions Tue Oct 2, 08:52

VnmrJ

File Edit View Experiments Help

2D Cryo ArrayedSpectra Protocols Frame Viewport 1D

Experiment Panel

Misc

PROTON	CARBON
(HH) gCOSY	(HC) gHMBCAD
zTOCSY	gHSQCAD
(C) DEPT	

Study Queue

Setup Show Time Acquire Stop Est. pw90 Sequence

Info/acquire proc/plot

Experiment: PROTON Obs: H1 Dec: C13 Sample: Eject Insert

Sample Info Use My Prefs

Sample Name: Alex_Bromination

Lot Number:

Notebook: Page:

Solvent: CDCI3

Comments: Position of bromination for Alex (PJP)

E-Address:

Current user: ijd Login Logout

Pre-Acquisition

Find z0

Grad. Shim

Acquisition

Scans: 16

Relax Delay: 1 sec

Spec width: -2.0 to: 14.0

Pulse Angle: 45 deg.

Post-Acquisition

Process

Plot

Display

Save

Tray auto_2007100

Display: LocQ

New Sample to Automation

StudyQ Action

Temp 300.0 C Spin 0 Hz Lock 0.1 Probe P8898_walkup

Idle

Locations accepted

[Terminal] vnmr500~vnmr1 :/hom VnmrJ

Vnmr500 Cheat Sheet

Remember: DayQ is 20 mins MAX
NightQ begins at 6pm

Any problems, see: Dr. Iain Day
2R113 / 2R115
tel: 6622

Click here to submit the next sample

The screenshot shows the VnmrJ software interface with the following components and callouts:

- Callout 1:** Points to the "Sample Name:" field in the "Sample Info" section.
- Callout 2:** Points to the "Solvent:" dropdown menu, which is currently set to "CDCl3".
- Callout 3:** Points to the "Comments:" text area.
- Callout 4:** Points to the "Next selection:" radio buttons, with "day" selected.
- Callout 5:** Points to the "Experiment Panel" on the left, specifically the "PROTON" experiment option.
- Callout 6:** Points to the "new Sample" entry in the "Study Queue" list.
- Callout 7:** Points to the grid of numbered circles (1-50) in the center, with circle 6 highlighted.
- Callout 8:** Points to the "Submit" and "Done" buttons in the "Add Next Selection to:" section.

Other visible interface elements include the "File Edit View Experiments Help" menu, the "2D Cryo ArrayedSpectra" tabs, the "Sample#: 3 (0 study)" status box, the "Experiment Queue" list showing "PROTON_01", and the bottom status bar with "Temp 30.0 C", "Spin 0 Hz", "Lock 0.1", and "Probe P8898_walkup".